

# Rapid sequencing DNA V14 - barcoding (SQK-RBK114.24 or SQK-RBK114.96)

Version: RBK\_9176\_v114\_revM\_27Nov2022  
Last update: 03/04/2024



Flow Cell Number: .....

DNA Samples: .....

## Before start checklist

### Materials

- ☐ 200 ng gDNA per sample
- ☐ Rapid Barcoding Kit 24 V14 (SQK-RBK114.24) OR Rapid Barcoding Kit 96 V14 (SQK-RBK114.96)

### Consumables

- ☐ MinION and GridION Flow Cell
- ☐ Qubit dsDNA HS Assay Kit (Invitrogen, Q32851)
- ☐ Bovine Serum Albumin (BSA) (50 mg/ml) (e.g. Invitrogen™ UltraPure™ BSA 50 mg/ml, AM2616)
- ☐ Freshly prepared 80% ethanol in nuclease-free water
- ☐ Nuclease-free water (e.g. ThermoFisher, AM9937)
- ☐ Qubit™ Assay Tubes (Invitrogen, Q32856)
- ☐ Eppendorf twin.tec® PCR plate 96 LoBind, semi-skirted (Eppendorf™, cat # 0030129504) with heat seals
- ☐ 2 ml Eppendorf DNA LoBind tubes
- ☐ 0.2 ml thin-walled PCR tubes
- ☐ 1.5 ml Eppendorf DNA LoBind tubes

### Equipment

- ☐ MinION or GridION device
- ☐ MinION and GridION Flow Cell Light Shield
- ☐ Microplate centrifuge, e.g. Fisherbrand™ Mini Plate Spinner Centrifuge (Fisher Scientific, 11766427)
- ☐ Timer
- ☐ Thermal cycler or heat blocks
- ☐ Magnetic rack
- ☐ Hula mixer (gentle rotator mixer)
- ☐ Ice bucket with ice
- ☐ Qubit fluorometer (or equivalent for QC check)
- ☐ Pipettes and pipette tips Multichannel, P2, P20, P100, P200, P1000

## INSTRUCTIONS

## NOTES/OBSERVATIONS

### Library preparation

Minimum Rapid Barcode use requirements

Check your flow cell.

- ☐ Program the thermal cycler: 30°C for 2 minutes, then 80°C for 2 minutes.
- Thaw kit components at RT, spin down briefly using a microfuge and mix by pipetting as indicated by the table below:
  - ☐ Rapid Barcodes (RB01-24) or Rapid Barcode Plate (RB01-96): not frozen, briefly spin down, mix well by pipetting
  - ☐ Rapid Adapter (RA): not frozen, briefly spin down, mix well by pipetting
  - ☐ AMPure XP Beads (AXP): thaw at RT, briefly spin down, mix by pipetting or vortexing immediately before use
  - ☐ Elution Buffer (EB): thaw at RT, briefly spin down, mix well by pipetting
  - ☐ Adapter Buffer (ADB): thaw at RT, briefly spin down, mix by vortexing

# Rapid sequencing DNA V14 - barcoding (SQK-RBK114.24 or SQK-RBK114.96)

Version: RBK\_9176\_v114\_revM\_27Nov2022  
Last update: 03/04/2024



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Prepare the DNA in Nuclease-free water.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Transfer 200 ng of genomic DNA per sample into 0.2 ml thin-walled PCR tubes or an Eppendorf twin.tec® PCR plate 96 LoBind.</li> <li><input type="checkbox"/> Adjust the volume of each sample to 10 µl with Nuclease-free water.</li> <li><input type="checkbox"/> Pipette mix the tubes for 10-15 times to avoid unwanted shearing</li> <li><input type="checkbox"/> Spin down briefly in a microfuge</li> </ul> <p>In the 0.2 ml thin-walled PCR tubes or an Eppendorf twin.tec® PCR plate 96 LoBind, mix the following:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 10 µl Template DNA (200 ng from previous step)</li> <li><input type="checkbox"/> 1.5 µl Rapid Barcodes (RB01-24 or RB01-96, one for each sample)</li> </ul> <ul style="list-style-type: none"> <li><input type="checkbox"/> Ensure the components are thoroughly mixed by pipetting and spin down briefly.</li> <li><input type="checkbox"/> Incubate the tubes or plate at 30°C for 2 minutes and then at 80°C for 2 minutes. Briefly put the tubes or plate on ice to cool.</li> <li><input type="checkbox"/> Spin down the tubes or plate to collect the liquid at the bottom.</li> <li><input type="checkbox"/> Pool all barcoded samples in a clean 2 ml Eppendorf DNA LoBind tube, noting the total volume.</li> <li><input type="checkbox"/> Resuspend the AMPure XP Beads (AXP) by vortexing.</li> </ul>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Ensure you have sufficient capacity in your reaction tube for all the reagents.</li> </ul>	
<p>Add an equal volume of resuspended AMPure XP Beads (AXP) to the entire pooled barcoded sample, and mix by flicking the tube.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Volume of AMPure XP Beads (AXP) added <ul style="list-style-type: none"> <li>- Volume per sample: 11.5 µl</li> <li>- For 4 samples: 46 µl</li> <li>- For 12 samples: 138 µl</li> <li>- For 24 samples: 276 µl</li> <li>- For 48 samples: 552 µl</li> <li>- For 96 samples: 1,000 µl</li> </ul> </li> <li><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 10 minutes at RT.</li> <li><input type="checkbox"/> Prepare at least 2 ml of fresh 80% ethanol in Nuclease-free water.</li> <li><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.</li> <li><input type="checkbox"/> Keep the tube on the magnet and wash the beads with 1 ml of freshly prepared 80% ethanol without disturbing the pellet. Remove the ethanol using a pipette and discard.</li> <li><input type="checkbox"/> Repeat the previous step.</li> <li><input type="checkbox"/> Briefly spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for 30 seconds, but do not dry the pellet to the point of cracking.</li> </ul>	

# Rapid sequencing DNA V14 - barcoding (SQK-RBK114.24 or SQK-RBK114.96)

Version: RBK\_9176\_v114\_revM\_27Nov2022  
Last update: 03/04/2024



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Remove the tube from the magnetic rack and resuspend the pellet in 15 µl Elution Buffer (EB) per 24 barcodes used.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Volume of Elution Buffer (EB) <ul style="list-style-type: none"> <li>- For 24 barcodes: 15 µl</li> <li>- For 48 barcodes: 30 µl</li> <li>- For 72 barcodes: 45 µl</li> <li>- For 96 barcodes: 60 µl</li> </ul> </li> <li><input type="checkbox"/> Incubate for 10 minutes at RT.</li> <li><input type="checkbox"/> Pellet the beads on a magnet until the eluate is clear and colourless, for at least 1 minute.</li> </ul> <p>Remove and retain the full volume of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Remove and retain the eluate which contains the DNA library in a clean 1.5 ml Eppendorf DNA LoBind tube</li> <li><input type="checkbox"/> Dispose of the pelleted beads</li> </ul>	
Quantify 1 µl of eluted sample using a Qubit fluorometer.	
<ul style="list-style-type: none"> <li><input type="checkbox"/> Transfer 11 µl of the sample into a clean 1.5 ml Eppendorf DNA LoBind tube.</li> </ul> <p>In a fresh 1.5 ml Eppendorf DNA LoBind tube, dilute the Rapid Adapter (RA) as follows and pipette mix:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 1.5 µl Rapid Adapter (RA)</li> <li><input type="checkbox"/> 3.5 µl Adapter Buffer (ADB)</li> <li><input type="checkbox"/> Add 1 µl of the diluted Rapid Adapter (RA) to the barcoded DNA.</li> <li><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</li> <li><input type="checkbox"/> Incubate the reaction for 5 minutes at RT.</li> </ul>	
The prepared library is used for loading into the flow cell. Store the library on ice until ready to load.	
<b>Priming and loading the MinION and GridION Flow Cell</b>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Please note, this kit is only compatible with R10.4.1 flow cells (FLO-MIN114).</li> </ul>	
<ul style="list-style-type: none"> <li><input type="checkbox"/> Thaw the Sequencing Buffer (SB), Library Beads (LIB) or Library Solution (LIS, if using), Flow Cell Tether (FCT) and Flow Cell Flush (FCF) at RT before mixing by vortexing. Then spin down and store on ice.</li> </ul>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> For optimal sequencing performance and improved output on MinION R10.4.1 flow cells (FLO-MIN114), we recommend adding Bovine Serum Albumin (BSA) to the flow cell priming mix at a final concentration of 0.2 mg/ml.</li> </ul>	
<p>To prepare the flow cell priming mix with BSA, combine the following reagents in a fresh 1.5 ml Eppendorf DNA LoBind tube. Mix by inverting the tube and pipette mix at RT:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 1,170 µl Flow Cell Flush (FCF)</li> <li><input type="checkbox"/> 5 µl Bovine Serum Albumin (BSA) at 50 mg/ml</li> <li><input type="checkbox"/> 30 µl Flow Cell Tether (FCT)</li> <li><input type="checkbox"/> 1,205 µl Final total volume in tube</li> </ul>	

# Rapid sequencing DNA V14 - barcoding (SQK-RBK114.24 or SQK-RBK114.96)

Version: RBK\_9176\_v114\_revM\_27Nov2022  
Last update: 03/04/2024



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<input type="checkbox"/> Open the MinION or GridION device lid and slide the flow cell under the clip. Press down firmly on the flow cell to ensure correct thermal and electrical contact.	
<input type="checkbox"/> Slide the flow cell priming port cover clockwise to open the priming port.	
<b>IMPORTANT</b> <input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 µl, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.	
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles:</p> <input type="checkbox"/> Set a P1000 pipette to 200 µl <input type="checkbox"/> Insert the tip into the priming port <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, to draw back 20-30 µl, or until you can see a small volume of buffer entering the pipette tip <p>Note: Visually check that there is continuous buffer from the priming port across the sensor array.</p> <input type="checkbox"/> Load 800 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for five minutes. During this time, prepare the library for loading by following the steps below. <input type="checkbox"/> Thoroughly mix the contents of the Library Beads (LIB) by pipetting.	
<b>IMPORTANT</b> <input type="checkbox"/> The Library Beads (LIB) tube contains a suspension of beads. These beads settle very quickly. It is vital that they are mixed immediately before use.	
<p>In a new 1.5 ml Eppendorf DNA LoBind tube, prepare the library for loading as follows:</p> <input type="checkbox"/> 37.5 µl Sequencing Buffer (SB) <input type="checkbox"/> 25.5 µl Library Beads (LIB) mixed immediately before use, or Library Solution (LIS), if using <input type="checkbox"/> 12 µl DNA library <p>Complete the flow cell priming:</p> <input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible. <input type="checkbox"/> Load 200 µl of the priming mix into the flow cell priming port (not the SpotON sample port), avoiding the introduction of air bubbles. <input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading. <input type="checkbox"/> Add 75 µl of the prepared library to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next. <input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port and close the priming port.	
<b>IMPORTANT</b> <input type="checkbox"/> Install the light shield on your flow cell as soon as library has been loaded for optimal sequencing output.	

Rapid sequencing DNA V14 - barcoding (SQK-RBK114.24 or SQK-RBK114.96)

Version: RBK\_9176\_v114\_revM\_27Nov2022  
Last update: 03/04/2024



Flow Cell Number: ..... DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Place the light shield onto the flow cell, as follows:</p> <ul style="list-style-type: none"><li><input type="checkbox"/> Carefully place the leading edge of the light shield against the clip. Note: Do not force the light shield underneath the clip.</li><li><input type="checkbox"/> Gently lower the light shield onto the flow cell. The light shield should sit around the SpotON cover, covering the entire top section of the flow cell.</li></ul>	
<p>Close the device lid and set up a sequencing run on MinKNOW.</p>	
<p><b>Flow cell reuse and returns</b></p>	
<ul style="list-style-type: none"><li><input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Flow Cell Wash Kit protocol and store the washed flow cell at 2-8°C.</li><li><input type="checkbox"/> Alternatively, follow the returns procedure to flush out the flow cell ready to send back to Oxford Nanopore.</li></ul>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"><li><input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol.</li></ul>	